

## AUDION FAMILY AMPLIFIER FEATURE MATRIX

FEATURES	Audion™ 4	Audion™ 6	Audion™ 8	Audion™ 16
Acous-TAP™				•
Adaptive Directionality		•	•	•
Audio Inputs, Number of	4	4	4	4
Compression Channels	4	6	8	16
Data Logging			•	•
Dual Biquad Filter				•
Dynamic Contrast Detection™	•	•	•	•
Feedback Managers:				
Adaptive Feedback Cancelling	•	•		
Reliant™ Adaptive Feedback Cancelling			3rd Generation	4th Generation
Fixed Directional Processing	•	•	•	•
Gain Adjustment Bands	12	12	12	16
High Power Compatible	•	•	•	•
In-situ Tone Generator		•	•	•
Input Filter Low Cut				•
Layered Noise Reduction™	10-Channel	12-Channel	12-Channel	12-Channel
Low Battery Warning Indicator	Programmable On/Off	Programmable On/Off	Programmable On/Off	Programmable On/Off
Low-level Expansion	•	•	•	•
M+T Input Mode	•	•	•	•
Maximum Power Output Limiting	Broadband	6-Channel	8-Channel	16-Channel
Program Switch Tones	Programmable On/Off	Programmable On/Off	Programmable On/Off	Programmable On/Off
Programmable Power-on Delay		•	•	•
Programmable Tone Generator		•	•	•
Programmable VC Range	•	•	•	•
Random Noise Generator			•	•
Tone Amplitude Adjustment	•	•	•	•
Tone Frequency Adjustment	•	•	•	•
Voice Prompts			•	•
Volume Control Type	Analog or Digital	Analog or Digital	Analog or Digital	Analog or Digital
Wide Dynamic Range Compression	•	•	•	•
Wind Noise Reduction				•
Hybrid Dimensions	.226" x .133"	.226" x .134"	.226" x .134"	.226" x .134"
Hybrid Mini Dimensions*				.210" x .116"

\*Mini Hybrids only have two inputs



Corporate Office:  
Germany Office:  
[www.IntriCon.com](http://www.IntriCon.com)  
© 2018 IntriCon Corp

1260 Red Fox Road, Arden Hills, Minnesota 55112  
Kesselschmiedstrasse 10, D-85354 Freising, Germany

23017\_B JUN-2018

TEL: 651-636-9770 FAX: 651-636-8944  
TEL: 49-8161-48040 FAX: 49-8161-480418

Contents Subject to Change Without Notice